

RESPIRONICS SIMPLYGO MINI

PORTABLE OXYGEN CONCENTRATOR

Preventive Maintenance

Periodic intervals:



Clean and replace the cannula and tubing as instructed by the manufacturer and your equipment provider.

When changing residents:



Clean and disinfect the exterior covers of the device and rechargeable battery between patient use and as needed.



#E1773 shown
with brown
backpack

PHILIPS
RESPIRONICS

- Delivers the highest oxygen output for portable concentrators weighing 5 lbs or less
- External battery design is easy to install and change
- LCD touchscreen is easy to read and use
- Two battery option available: standard (up to 4.5 hrs) and extended (up to 9 hrs)

Specifications:	
Flow Rate:	0.5 - 3.0 lpm in 0.5 lpm increments
O ₂ Purity %:	90% +5.5/-3%
Sound Level:	43dBA typical at pulse setting 2"
Battery Life:	4½ hours
Dimensions:	83/10"W x 33/5"D x 92/5"H
Weight:	5 lbs.
Warranty:	3-year warranty on the system, 1-year warranty on the sieve