



**Keep your
wheelchairs rolling
smoothly for years**

We have the quality parts you need to help keep your chairs in great shape! No matter what part needs to be replaced, you're sure to find the replacement within our thousands of options.

Call your account manager anytime for assistance.

Panacea® Heavy-Duty Reclining Wheelchair

Replacement Parts

Upholstery

- 1a. **#B0121** Headrest Pillow, 22"
- 1b. **#B0122** Headrest, 22", without Pillow
- 1c. **#B0123** Headrest, Upholstery Only, 22"
- 1d. **#B0124** Headrest Post
- 2. **#B0126** Back Upholstery, 22"
- 2a. **#A2623** Back Upholstery Screws
- 17. **#B0127** Seat Upholstery, 22"
- 19. **#A2624** Seat Upholstery Screws
- 30. **#B0125** Backrest Cane
- 34. **#B0128** Seat Saddle

Frame

- 23. **#B0146** Slide Tube
- 28. **#B0143** Cross Frame, 22"
- 28. **#A2635** Cross Frame Bolt (not pictured)
- 29. **#B0144** Side Frame (Rt)
- 29. **#B0145** Side Frame (Lt)

Footrest

- 25. **#B0153** Elevating Legrests
- 25a. **#B0152** Calf Pad Black Vinyl
- 26. **#B0151** Footplate
- 27. **#A1842** Footrest Bumper

Armrest

- 18. **#A1802** Seat Rail Guide Detach Arm, Rear
- 20. **#B0147** Desk Arm Assembly (Lt)
- 20. **#B0148** Desk Arm Assembly (Rt)
- 20. **#B0563** Adj. Height Desk Arm Assembly (Lt) (not pictured)
- 20. **#B0562** Adj. Height Desk Arm Assembly (Rt) (not pictured)
- 20. **#B0149** Full Arm Assembly (Lt) (not pictured)
- 20. **#B0150** Full Arm Assembly (Rt) (not pictured)
- 20. **#B0565** Adj. Height Full Arm Assembly (Lt) (not pictured)
- 20. **#B0564** Adj. Height Full Arm Assembly (Rt) (not pictured)
- 22. **#A1799** Desk Arm Pad
- 22. **#A1800** Full Arm Pad (not pictured)
- 24. **#A2641** Arm Release Button

Reclining Functionality

- 21. **#A2638** Handgrip
- 31. **#B0154** Recline Release Lever (Rt) (No Cable)
- 31. **#B0155** Recline Release Lever (Lt) (No Cable)
- 32. **#A1875** Cable
- 33. **#B0156** Reclining Cylinder (Rt)
- 33. **#B0157** Reclining Cylinder (Lt)

Wheels

- 3. **#B0129** Rear Wheel/Chrome Handrim
- 4. **#B0132** Rear Wheel Bearing
- 5. **#B0131** Rear Axle
- 6. **#B0130** Hub Cap
- 7. **#B0133** Anti-Tippers without Wheels
- 8. **#A2628** Chrome Handrim
- 9. **#B0134** Handrim Screws & Spacers
- 10. **#B0137** Brake Tip
- 11. **#B0135** Brake (Rt)
- 11. **#B0136** Brake (Lt)
- 12. **#B0141** Front Axle
- 13. **#B0140** Caster Bearing
- 14. **#B0139** Fork
- 15. **#B0142** Fork Stem Bearing
- 16. **#A1840** Dust Cap
- 35. **#B0138** Caster