

# DRIVE DEVILBISS® IGO2

## PORTABLE OXYGEN CONCENTRATOR



#347YR,  
#347YT



- SmartDose Auto-Adjusting Conserver Technology included
- Settings 1-5 pulse dose
- Rugged design, including protective overmold
- Weighs less than 5 lbs
- Most sensitive conserver trigger on the POC market

## Preventive Maintenance

### Periodic intervals:



Clean and replace the cannula and tubing as instructed by the manufacturer and your equipment provider.

### When changing residents:



Clean and disinfect the exterior covers of the device and rechargeable battery between patient use and as needed.

### Specifications:

Flow Rate:	1 - 5 pulse settings
O <sub>2</sub> Purity %:	90% +4/-3%
Sound Level:	37 dBA at setting 2
Battery Life:	3.5 hours on pulse setting 2
Dimensions:	3.5"W x 8.6"D x 8.4"H
Weight:	4.95 lbs.
Warranty:	5 years device, 3 years compressor, 1 year sieve beds and battery