

# AIRSEP VISIONAIRE

## 5-LITER OXYGEN CONCENTRATOR

### Preventive Maintenance

#### Once Every 3 Years:



Verify the oxygen purity, flow rate and pressure of your concentrator by using an oxygen analyzer.

#### When changing patients:



Discard oxygen tubing, cannula/mask, oxygen outlet connector, and humidifier bottle. Replace cabin filter. Clean concentrator cabinet. Inspect plugs, cords and components. Check oxygen concentration.



#A9970



A Chart Industries Company

- No filters to change; easiest model to maintain
- Optional alarms for purity and flow rate
- Economical price point

#### Specifications:

Flow Rate:	1 - 5.0 lpm
O <sub>2</sub> Purity %:	87% - 95%
Power Consumption:	290W
Sound Level:	40 dB
Dimensions:	14 <sup>1</sup> / <sub>16</sub> "W x 11 <sup>1</sup> / <sub>2</sub> "D x 20 <sup>1</sup> / <sub>8</sub> "H
Weight:	30 lbs.
Warranty:	3 years*

\*2-year extended warranty available (#C5939)