



Replacement Parts for
Panacea®
Reclining Wheelchair

Keep your wheelchairs rolling smoothly for years

We have the quality parts you need to help keep your chairs in great shape! No matter what part needs to be replaced, you're sure to find the replacement within our thousands of options.

Call your account manager anytime for assistance.

Panacea® Reclining Wheelchair

Replacement Parts

Armrest/Footrest/Headrest

- 1. **#A1859** Headrest Pillow
- 2. **#A1856** Headrest, 16"
- 2. **#A1857** Headrest, 18"
- 2. **#A1858** Headrest, 20"
- 20. **#A1868** Elevated Legrest, Pair
- 26. **#A2651** Desk Arm
- 26a. **#A1799** Desk Arm Pad
- 26a. **#A1800** Full Arm Pad

Upholstery

- 3. **#A1853** Back Upholstery, 16"
- 3. **#A1854** Back Upholstery, 18"
- 3. **#A1855** Back Upholstery, 20"
- 24. **#A1850** Seat Upholstery, 16"
- 24. **#A1851** Seat Upholstery, 18"
- 24. **#A1852** Seat Upholstery, 20"
- 25. **#A2580** Upholstery Screw for Back

Frame

- 4. **#A1872** Backrest Cane
- 5. **#A2650** Handgrip
- 6. **#A1878** Release Handle/Cable, Right
- 6. **#A1879** Release Handle/Cable, Left
- 7. **#A1873** Spreader Bar Knob
- 8. **#A2645** Spreader Bar, 16"
- 8. **#A2646** Spreader Bar, 18"
- 8. **#A2647** Spreader Bar, 20"
- 9. **#A1875** Cable
- 10. **#A1876** Reclining Piston, Right
- 10. **#A1877** Reclining Piston, Left
- 14. **#A1867** Anti-Tippers
- 27. **#A2649** Side Frame
- 28. **#A3600** Cross Frame, 16"
- 28. **#A3601** Cross Frame, 18"
- 28. **#A3602** Cross Frame, 20"
- 29. **#A1874** Flip Lock-Arm Release
- 30. **#A1871** Slide Tube

Wheels

- 11. **#A1788** Handrim Screws & Spacers, Set of 9
- 12. **#A2643** Rear Axle Bolt
- 12a. **#A1862** Hub Cap
- 13. **#A1863** Rear Bearing, Set of 4
- 15. **#A1865** Rear Wheel
- 23. **#A1860** Wheel Lock, Right
- 23. **#A1861** Wheel Lock, Left

Casters

- 16. **#A1866** Caster Bearing
- 17. **#A2644** Caster Axle
- 18. **#A1864** Caster
- 19. **#A2648** Fork
- 21. **#A1869** Fork Stem Bearings
- 22. **#A1870** Dust Cap