

START HERE

DS smart™ Hardware Kit Contents

Thank you for your purchase of the DS smart™ platform for use with the Attendant® Connected Vital Signs Monitor. Included in this package are all the materials you'll need to get started.

Box Contents:

- Tablet, tablet case, tablet mount, power strip assembly, cable jacket, hex key, USB cable and stylus tether
- **HARDWARE ASSEMBLY:** This is the first step for assembly and installation of DS smart™. Take contents of box contents and either:
 - Watch assembly video located at <http://bit.ly/TabletMountAssembly> OR
 - Use enclosed **HARDWARE ASSEMBLY GUIDE** document for instructions. Note you will need the enclosed hex key, a small flathead screwdriver and a medium Phillips screwdriver to assemble the hardware.
- **TABLET SETUP:** Once the hardware has been assembled, take the tablet (in the tablet case) and use the enclosed **INITIAL TABLET SETUP INSTRUCTIONS** document to install the DS smart™ app.
- Troubleshooting Guide (includes alligator clip to attach to basket on vital signs monitor cart)

